

TransScript® Green miRNA Two-Step qRT-PCR SuperMix

Cat.No. AQ202

Storage: at -20°C for two years

Description

TransScript® Green miRNA Two-Step qRT-PCR SuperMix provides all the necessary components for the detection and quantification of miRNA from total RNA, small RNA or other miRNA containing samples. TransScript® miRNA RT Enzyme Mix (Poly(A) polymerase and reverse transcriptase are included) and 2×TS miRNA Reaction Mix are supplied to efficiently add Poly(A) tails and synthesize first-strand cDNA. PerfectStart™ Green qPCR SuperMix is used for miRNA quantification.

Highlights

- Optimal ratio of Poly(A) polymerase to reverse transcriptase and optimized reaction buffer, so as to ensure high transcription efficiency of miRNA.
- One-step completion of Poly(A) tail addition and first-strand cDNA synthesis in a single tube.
- High amplification efficiency, specificity and sensitivity ensured by PerfectStart™ Green qPCR SuperMix, leading to accurate data.
- Passive reference dyes compatible with different qPCR instruments (normalize the fluorescent signals between reactions).

Application

Multiple copy and low copy gene detection

Passive Reference Dye

- Passive Reference Dye I (50×)
ABI Prism® 7000/7300/7700/7900, ABI Step One®, ABI Step One Plus®
- Passive Reference Dye II (50×)
ABI Prism® 7500, ABI Prism® 7500 Fast, ABI Q6, ABI QuantStudio® 6/7 Flex, ABI ViiA® 7, Stratagene Mx3000®/Mx3005P®, Qiagen Corbett Rotor-Gene® 3000
- No Passive Reference Dye
Roche LightCycler® 480, Roche Light Cycler® 96, MJ Research Chromo4®, MJ Research Opticon® 2, Takara TP-800®, Bio-Rad iCycler iQ®, Bio-Rad iCycler iQ5®, Bio-Rad CFX96®, Bio-Rad C1000® Thermal Cycler, Thermo Scientific Pikoreal® 96, Qiagen Corbett Rotor- Gene® 6000, Qiagen Corbett Rotor-Gene® G, Qiagen Corbett Rotor-Gene® Q

Kit Contents

Component	AQ202-01
TransScript® miRNA RT Enzyme Mix	20 µl
2×TS miRNA Reaction Mix	200 µl
Universal miRNA qPCR Primer (10 µM)	200 µl
2×PerfectStart™ Green qPCR SuperMix	5×1 ml
Passive Reference Dye (50×)	200 µl
RNase-free Water	1 ml

Tail addition and first-strand cDNA synthesis

1. Reaction Components

Component	Volume
Total RNA/ miRNA*	x μ l
<i>TransScript</i> [®] miRNA RT Enzyme Mix	1 μ l
2 \times TS miRNA Reaction Mix	10 μ l
RNase-free Water	to 20 μ l

* Total RNA \leq 5 μ g. Since miRNA cannot be directly quantified by spectrophotometer, we suggest using 1-9 μ l for 20 μ l reaction.

2. Mix gently, and incubate at 37°C for 1 hour.

3. Incubate at 85°C for 5 seconds to inactivate RT Enzyme Mix.

Suggested qPCR conditions (20 μ l reaction volume)

Component	Volume	Final Concentration
cDNA* ¹	Variable	as required
Forward Primer (10 μ M)* ²	0.4 μ l	0.2 μ M
Universal miRNA qPCR Primer (10 μ M)	0.4 μ l	0.2 μ M
2 \times <i>PerfectStart</i> [™] Green qPCR SuperMix	10 μ l	1 \times
Passive Reference Dye (50 \times) (optional)	0.4 μ l	1 \times
Nuclease-free Water	Variable	-
Total Volume	20 μ l	-

*1. We suggest diluting the synthesized cDNA 5-10 folds.

*2. Upstream primer is target miRNA specific primer, which will be designed by customers according to target miRNA.

Thermal cycling conditions (three-step)

94°C 30 sec
 94 °C 5 sec
 50-60°C 15 sec★
 72°C 10 sec★

40-45 cycles

Dissociation Stage

Thermal cycling conditions (two-step)

94°C 30 sec
 94°C 5 sec
 60°C 30 sec★

40-45 cycles

Dissociation Stage

Fluorescent signals can be collected during the annealing or extension stage. For ABI qPCR instrument, we suggest using the following signal collecting time:

* For ABI Prism[®] 7700/7900, the time to 30 seconds.

* For ABI Prism[®] 7000/7300, the time to 31 seconds.

* For ABI Prism[®] 7500, the time to 34 seconds.

* For ABI ViiA[®] 7, the time is at least 19 seconds.

Two-step qPCR is more suitable for higher specificity assay.

Three-step qPCR is more suitable for higher sensitivity assay.

FOR RESEARCH USE ONLY