

TransStart® Green qPCR SuperMix UDG

Cat. No. AQ111

Storage: at -20°C in dark for two years

Description

TransStart® Green qPCR SuperMix UDG is a ready-to-use qPCR cocktail containing all components, except primer and template. It contains TransStart® Taq DNA Polymerase, UDG, SYBR Green I, dNTPs, PCR enhancer and stabilizer. qPCR SuperMix is provided at 2× concentration and can be used at 1× concentration by adding template, primer, passive reference dye (optional) and Nuclease-free Water.

Highlights

- TransStart® Taq DNA Polymerase, hot start with double blocking technique, improves sensitivity, enhances specificity and generates more accurate data.
- Double cation (K^+ , NH_4^+) buffer enhances specificity and reduces primer-dimer formation.
- Passive reference dyes are provided for different qPCR instruments.
- UDG and dUTP avoid cross contamination.

Passive Reference Dye

- Passive Reference Dye I (50×)

ABI Prism® 7000/7300/7700/7900, ABI Step One®, ABI Step One Plus®

- Passive Reference Dye II (50×)

ABI Prism® 7500, ABI Prism® 7500 Fast, ABI Q6, ABI QuantStudio® 6/7 Flex, ABI ViiA® 7, Stratagene Mx3000® /Mx3005P®, Qiagen Corbett Rotor-Gene® 3000

- No Passive Reference Dye

Roche LightCycler® 480, Roche Light Cycler® 96, MJ Research Chromo4®, MJ Research Opticon® 2, Takara TP-800®, Bio-Rad iCycler iQ®, Bio-Rad iCycler iQ5®, Bio-Rad CFX96®, Bio-Rad C1000® Thermal Cycler, Thermo Scientific Pikoreal® 96, Qiagen Corbett Rotor- Gene® 6000, Qiagen Corbett Rotor-Gene® G, Qiagen Corbett Rotor-Gene® Q

Kit Contents

Component	AQ111-01	AQ111-02	AQ111-03
TransStart® Green qPCR SuperMix UDG (2×)	1 ml	5×1 ml	15×1 ml
Passive Reference Dye (50×)	40 µl	200 µl	600 µl
Nuclease-free Water	1 ml	5 ml	3×5 ml

Reaction Components (20 µl)

Component	Volume	Final Concentration
Template	Variable	as required
Forward Primer (10 µM)	0.4 µl	0.2 µM
Reverse Primer (10 µM)	0.4 µl	0.2 µM
2×TransStart® Green qPCR SuperMix UDG	10 µl	1×
Passive Reference Dye (50×) (optional)	0.4 µl	1×
Nuclease-free Water	Variable	-
Total Volume	20 µl	-

For genomic DNA, we suggest using 1 pg-1 µg template; for plasmid DNA, we suggest using 10-10⁷ copies.



Thermal cycling conditions (three-step)

50°C 2 min (UDG Incubation)
94°C 10 min (UDG Inactivation)
94°C 5 sec
50-60°C 15 sec*
72°C 10 sec* } 40-45 cycles
Dissociation Stage

Thermal cycling conditions (two-step)

50°C 2 min (UDG Incubation)
94°C 10 min (UDG Inactivation)
94°C 5 sec
60°C 30 sec* } 40-45 cycles
Dissociation Stage

Fluorescent signals can be collected during the annealing or extension stage. For ABI qPCR instrument, we suggest using the following signal collecting time:

- * For ABI Prism® 7700/7900, the time to 30 seconds.
- * For ABI Prism® 7000/7300, the time to 31 seconds.
- * For ABI Prism® 7500, the time to 34 seconds.
- * For ABI ViiA® 7, the time is at least 19 seconds.

Two-step qPCR is more suitable for higher specificity assay.

Three-step qPCR is more suitable for higher sensitivity assay.

Note

Completely thaw the contents in the tube and mix well before each use.

For research use only, not for clinical diagnosis.

Service telephone +86-10-57815020

Service email complaints@transgen.com

