

EasyPure® FFPE Tissue Genomic DNA Kit

Please read the manual carefully before use

Cat. No. EE191

Storage: At room temperature (15-25°C) for one year.

Description

EasyPure® FFPE Tissue Genomic DNA Kit is designed for extracting genomic DNA from formalin-fixed, paraffin-embedded (FFPE) tissue blocks and sections. It uses xylene deparaffinization method to remove paraffin, followed by enzymatic tissue lysis and a specialized incubation step to effectively reverse formaldehyde-induced cross-linking. DNA is specifically adsorbed using a silica-based column. The purified DNA is suitable for downstream applications including PCR, qPCR.

Features

- Resistant to formaldehyde inhibition, with strong lysis capability, fast extraction speed, and high yield.
- The spin column efficiently and specifically adsorbs DNA, effectively removing impurities such as proteins and salts, resulting in high-purity DNA.

Kit Contents

| Components | EE191-01 (50 rxns) |
|---|--------------------|
| Lysis Buffer 15 (LB15) | 11 ml |
| Binding Buffer 15 (BB15) | 11 ml |
| Clean Buffer 15 (CB15) | 6 ml |
| Wash Buffer 15 (WB15) | 12 ml |
| Elution Buffer (EB) | 25 ml |
| Proteinase K (20 mg/ml) | 1 ml |
| Genomic Spin Column with Collection Tubes | 50 each |

Before starting

- Add 100% ethanol to CB15 and WB15

| | Volume of 100% ethanol to add |
|------------------------|-------------------------------|
| Clean Buffer 15 (CB15) | 24 ml |
| Wash Buffer 15 (WB15) | 48 ml |

Procedures (all centrifugation steps are performed at room temperature)

1. Preparing materials

A. Tissue fixed in formalin, etc.: Place 10-30 mg chopped tissue into a 1.5 ml microcentrifuge tube.

- Add 500 µl 1×PBS and mix thoroughly by vortexing. Incubate at room temperature for 1 minute. Centrifuge at 12,000×g for 30 seconds and remove the supernatant.
- Repeat step a twice.

B. Paraffin-embedded tissue block: Trim excess paraffin off the tissue and scrap 10-30 mg tissue by a scalpel.

Paraffin-embedded tissue sections: Scrap tissue form 3-10 sections (5-10 µm thick) and place into a 1.5 ml microcentrifuge tube.

- Add 1 ml xylene to the sample, close the lid and vortex vigorously for 10 seconds. Centrifuge at 12,000×g for 2 minutes, then remove the supernatant by pipetting.

It is recommended to carefully perform this step in a fume hood because xylene is highly toxic chemical. Avoid contact with skin, eyes, and respiratory tract. Additionally, exercise caution to keep away from open flames during operation.

- Add 1 ml absolute ethanol to the pellet, and mix by vortexing. Centrifuge at 12,000×g for 2 minutes, then remove the supernatant by pipetting.
- Open the tube and incubate at room temperature or up to 37°C until all residual ethanol has evaporated.



2. Add 200 μ l of LB15 and 20 μ l of Proteinase K to resuspend the pellet, mix thoroughly by vortexing. Incubate at 56°C for 1 hour (or until the sample has been completely lysed).
3. Incubate at 90°C for 1 hour. Briefly centrifuge the tube to remove drops from the lid.
The incubation at 90°C can partially reverse formalin-induced crosslinking of nucleic acids. Longer incubation time and higher temperature may result in more fragmented DNA. Therefore, the heating device must be preheated to 90°C before incubation begins.
Optional: If RNA-free genomic DNA is needed, add 10 μ l of RNase A (20 mg/ml, Cat. No. GE101) to the sample after it has cooled to room temperature, incubate at room temperature for 2 minutes.
4. Add 200 μ l of BB15 to the sample, and mix thoroughly by vortexing (white precipitates may form).
5. Add 250 μ l of absolute ethanol, and mix thoroughly by vortexing (if there are white precipitates, vortex until they disappear). Briefly centrifuge the 1.5 ml tube and transfer all the lysate to a Genomic Spin Column, then centrifuge the tube at 12,000 \times g for 1 minute, discard the flow-through.
6. Add 500 μ l of CB15 (make sure that ethanol has been absolute ethanol), and centrifuge at 12,000 \times g for 30 seconds, discard the flow-through.
7. Add 500 μ l of WB15 (make sure that ethanol has been absolute ethanol), and centrifuge at 12,000 \times g for 30 seconds, discard the flow-through.
8. Repeat step 7 once.
9. Centrifuge the Spin Column at 15,000 \times g for 2 minutes. Open the tube and air-dry the Spin Column at room temperature to volatilize the residual ethanol.
10. Place the Spin Column in a clean 1.5 ml microcentrifuge tube. Add 30-100 μ l of Elution Buffer (for higher yield, prewarm the buffer to 65°C) or distilled water (pH >7.0) to the center of the membrane. Incubate at room temperature for 1 minute. Centrifuge at 12,000 \times g for 1 minute to elute the DNA.

Notes

- Too many starting materials will affect DNA yield.
- To ensure the quality of the extracted DNA, it is critical that the tissue is fresh prior to fixation. Use a fixation time of 14-24 hours (longer fixation times lead to more fragmented DNA). Thoroughly dehydrate samples prior to embedding because residual formalin can inhibit the Proteinase K digestion.
- Use sterile centrifuge tubes and pipette tips to avoid DNase contamination.
- Prolonged sample storage time (>2 years) may lead to extraction failure.
- It is recommended to measure the DNA concentration prior to PCR. Using 20 ng of DNA as template (in a 20 μ l reaction system) typically yields optimal amplification.

For research use only, not for clinical diagnosis.

Service telephone +86-10-57815020

Service email custserv@transgenbiotech.com

