

## *pEASY*<sup>®</sup>-Blunt Cloning Kit

Please read the user manual carefully before use.

**Cat. No.** CB101

### Storage

*Trans*1-T1 Phage Resistant Chemically Competent Cell at -70°C or below for six months; others at -18°C or below for nine months

### Description

*pEASY*<sup>®</sup>-Blunt Cloning Vector containing the *LacZ* gene allows for blue-white screening on plates supplemented with IPTG and X-gal. It is suitable for blunt-end cloning.

### Features

- Fast: 5 minutes.
- Simple: Just add the fragment.
- High cloning efficiency: High positive cloning rate.
- Kanamycin and Ampicillin resistance genes for selection.
- M13 Forward Primer, M13 Reverse Primer and T7 Promoter Primer for sequencing.
- T7 promoter for *in vitro* transcription.
- *Trans*1-T1 competent cells have high transformation efficiency and rapid growth rate, ensuring the number of clones and saving screening time.

### Kit Contents

Component	CB101-01 (20 rxns)	CB101-02 (60 rxns)
<i>pEASY</i> <sup>®</sup> - Blunt Cloning Vector (10 ng/μl)	20 μl	3×20 μl
Control Template (5 ng/μl)	5 μl	5 μl
Control Primers (10 μM)	5 μl	5 μl
M13 Forward Primer (10 μM)	50 μl	150 μl
M13 Reverse Primer (10 μM)	50 μl	150 μl
<i>Trans</i> 1-T1 Phage Resistant Chemically Competent Cell	10×100 μl	30×100 μl

### Preparation of PCR Products

1. Primer requirement: primer cannot be phosphorylated
2. PCR Enzyme: High-fidelity DNA polymerases that generate blunt-end products, such as *FastPfu* or *KD Plus* DNA Polymerase
3. Reaction conditions: in order to ensure the integrity of amplification products, 5-10 minutes of post-extension step is required.

After amplification, use agarose gel electrophoresis to verify the quality and quantity of PCR product

### Setting Up the Cloning Reaction System

Add following components into a microcentrifuge tube.

Component	Volume
PCR Product	0.5-4 μl
<i>pEASY</i> <sup>®</sup> -Blunt Cloning Vector	1 μl

Gently mix and incubate the mixture at room temperature (20°C-37°C) for 5 minutes, and then place the tube on the ice.

1. Optimal amount of insert

Molar ratio of vector to insert = 1:7 (1 kb, ~20 ng; 2 kb, ~40 ng)



2. Optimal volume of vector: 1  $\mu$ l

3. Optimal reaction volume: 3~5  $\mu$ l

If the volume is insufficient, add sterile water to adjust the final volume.

4. Optimal incubation time

(1) 0.1~1 kb (including 1 kb): 5~10 minutes\*

(2) 1~2 kb (including 2 kb): 10~15 minutes\*

(3) 2~3 kb (including 3 kb): 15~20 minutes\*

(4) For fragments larger than 3 kb: 20~30 minutes.

\*Use the maximum incubation time if the insert is purified using gel extraction.

• The optimal incubation temperature is 25°C. For fragment with high GC content, the reaction can be performed at 37°C. (It is recommended to use a PCR instrument for temperature control.)

### Transformation

1. Add the ligated products to 50  $\mu$ l of *Trans1*-T1 Phage Resistant Chemically Competent Cell (Cells are used immediately after thawing) and mix by tapping the tube gently.

2. Incubate on ice for 20~30 minutes.

3. Heat-shock the cells at 42°C for 30 seconds.

4. Immediately place the tube on ice for 2 minutes.

5. Add 250  $\mu$ l of room temperature SOC or LB medium. Shake the tube at 37°C (200 rpm) for 1 hour.

6. In the meantime, mix 8  $\mu$ l of 500 mM IPTG with 40  $\mu$ l of 20 mg/ml X-gal. Spread them evenly onto a selective LB plate. Place the plate at 37°C for 30 minutes.

7. Spread 200  $\mu$ l of cells on plate and incubate overnight. (If more colonies are desired, centrifuge at 1,500 $\times$ g for 1 minute, discard some of the supernatant, retain 100-150  $\mu$ l, resuspend the bacteria by gently tapping the tube, spread the entire bacterial liquid onto the plate, and incubate overnight.)

### Identification of Positive Clones and Sequencing

#### Analysis of positive clones

1. Pick a white single colony into 10  $\mu$ l of sterile water and vortex to mix.

2. Use 1  $\mu$ l of the mixture as template for 25  $\mu$ l PCR using M13 forward and M13 reverse primers.

3. PCR reaction conditions

94°C 10 min

94°C 30 sec

55°C 30 sec

72°C x min\*

72°C 5-10 min

30 cycles

\* (depends on the insert size and PCR enzymes) the PCR product size from vector self-ligation is 200 bp.

4. Analyze positive clones by restriction enzyme digestion and DNA sequencing.

Inoculate positive clones on LB/Amp<sup>+</sup> or LB/Kan<sup>+</sup> liquid medium, grow at 37°C for 6 hours at 200 rpm. Isolate plasmid DNA by plasmid MiniPrep Kit. Analyze plasmids by restriction enzyme digestion with proper restriction endonuclease.

### Sequencing

Analyze the sequence by sequencing with M13 F, M13 R and T7 promoter.

#### PCR for control insert (700 bp)

Component	Volume	Final Concentration
Control Template	1 $\mu$ l	0.1 ng/ $\mu$ l
Control Primers (10 $\mu$ M)	1 $\mu$ l	0.2 $\mu$ M
2 $\times$ EasyPfu PCR SuperMix	25 $\mu$ l	1 $\times$
Nuclease-free Water	Variable	-
Total Volume	50 $\mu$ l	-

PCR

94°C 2-5 min

94°C 30 sec


55°C 30 sec

72°C 1 min

72°C 10 min

30 cycles





**For research use only, not for clinical diagnosis.**

Service telephone +86-10-57815020

Service email [custserv@transgenbiotech.com](mailto:custserv@transgenbiotech.com)

