

## MagicPure® mRNA Kit

Please read the manual carefully before use.

Cat. No. EC511

Version No. Version 6.0

Storage: at 2-8°C for one year

### Description

MagicPure® mRNA Kit uses oligo(dT)-conjugated magnetic beads to specifically bind to poly(A) tailed mRNA. It is suitable for isolating mRNA from purified highly intact total RNA (0.01-10 µg, RIN value≥8). The isolated mRNA can be used in RT-PCR, qRT-PCR, next generation sequencing, or other applications. This kit is compatible with magnetic-rod high-throughput nucleic acid extractor.

### Highlights

- High-yield and high-purity isolated mRNA
- Simple workflow

### Sample Requirements

0.01-10 µg of purified highly intact total RNA (RIN value≥8)

### Kit Contents

Component	EC511-01 (24 rxns)	EC511-02 (96 rxns)
Binding Buffer 33 IV (BB33 IV)	1.3 ml	5 ml
Clean Buffer 33 (CB33)	1.3 ml	5 ml
Wash Buffer 33 III (WB33 III)	10 ml	40 ml
RNase-free Water	1.3 ml	5 ml
mRNA Beads IV	1.3 ml	5 ml

### Procedures

1. Take all reagents out from 2-8°C and **equilibrate to room temperature for 30 minutes. Mix well by vortexing.**
2. Prepare RNA sample: dilute the total RNA to 50 µl with RNase-free Water in a centrifuge tube.
3. Pipet 50 µl of mRNA Beads IV to the RNA sample. Mix well by pipetting up and down.
4. Place the centrifuge tube in a PCR machine at 65°C for 5 minutes, 25°C for 5 minutes, and hold at 4°C to allow the mRNA to bind to the magnetic beads.  
Note: Make sure that the beads have been mixed thoroughly prior to reaction.
5. Place the centrifuge tube on a magnetic stand for 5 minutes. Discard the supernatant carefully and completely.
6. Remove the centrifuge tube from the magnetic stand. Add 200 µl of WB33 III and mix well by pipetting. Place on the magnetic stand for 5 minutes and discard the supernatant carefully and completely.
7. Remove the centrifuge tube from the magnetic stand. Add 50 µl of CB33 and resuspend the beads by pipetting.
8. Heat the tube for 2 minutes at 80°C and cool to 25°C.
9. Add 50 µl of BB33 IV and mix well by pipetting. Incubate at room temperature for 5 minutes.
10. Place the centrifuge tube on the magnetic stand for 5 minutes. Discard the supernatant carefully and completely.
11. Remove the centrifuge tube from the magnetic stand. Add 200 µl of WB33 III and mix well by pipetting. Place on the magnetic stand for 5 minutes and discard the supernatant carefully and completely.
12. Choose the processing method according to the experimental process:
  - a. Use automated instrument is used for purification or the purified product is used for reverse transcription: Remove the tube from the magnetic stand. Add 18 µl RNase-free Water. Pipet 6 times to mix thoroughly. Incubate at 80°C for 2 minutes. Immediately place the tube on the magnetic stand for 5 minutes. After the solution is clear, carefully pipet 16 µl of the supernatant to into a new RNase-free centrifuge tube.
  - b. Purified product is used for RNA library construction, such as TransNGS Fast RNA Seq Library Prep Kit for Illumina Rapid RNA Library Prep Kit (Cat. No. KP701): Add 1×RNA Fragmentation Buffer according to the manual for library construction.
13. The purified product can be placed on ice for NGS library construction or other analysis applications (it is recommended to carry out subsequent reactions immediately), or it can be stored at -80°C.



#### Notes

- Please use RNase-free centrifuge tubes;
- The total RNA sample should be highly intact (RIN value>8). Otherwise the mRNA information will be partially lost.

**For research use only, not for clinical diagnosis**

Version number: V6.0-202505

Service telephone +86-10-57815020

Service email [complaints@transgen.com](mailto:complaints@transgen.com)

