

TransNGS® Library Amplification SuperMix

Please read the manual carefully before use.

Cat. No. KA101

Version No. Version 5.0

Storage: at -20°C for two years

Description

This product contains *TransStart® FastPfu* Fly DNA Polymerase, dNTPs and a reaction buffer optimized for library amplification. The concentration is 2×, which has the characteristics of high sensitivity, low preference and high fidelity. This product is suitable for the amplification of next-generation sequencing libraries. Under the recommended amplification conditions, it can achieve high fidelity and low preference during the amplification of next-generation sequencing libraries, and can meet the amplification of libraries with different GC contents. During amplification, just add template, primers and water to make the concentration of SuperMix solution 1×, and the amplified product is blunt-ended.

Highlights

- High fidelity amplification.
- Low amplification bias.
- High sensitivity and high specificity.
- Hot start.

Applications

- Next-generation sequencing library amplification.

Kit Contents

| Component | KA101-01 | KA101-02 |
|--|----------|----------|
| TransNGS® Library Amplification SuperMix | 1 ml | 5×1 ml |
| Nuclease-free Water | 1 ml | 5 ml |

Reaction Components

| Component | Volume | Final Concentration |
|--|----------|---------------------|
| Adapter-ligated DNA | Variable | Variable |
| Library Amplification Forward Primer (10 μM) | 2.5 μl | 0.5 μM |
| Library Amplification Reverse Primer (10 μM) | 2.5 μl | 0.5 μM |
| TransNGS® Library Amplification SuperMix | 25 μl | 1× |
| Nuclease-free Water | Variable | - |
| Total volume | 50 μl | - |

Recommended thermal cycling conditions

| | | |
|--------|--------|-----------------|
| 98°C | 3 min | } 2-15 cycles** |
| 98°C | 30 sec | |
| x°C* | 30 sec | |
| 72°C | 30 sec | |
| 72°C | 3 min | |
| ≤ 10°C | Hold | |

*Depending on the PCR primer length and GC contents.

**Refer to the manual of library construction kit for the number of cycles.



Notes

- All components should be thawed and mixed thoroughly before use.
- We suggest to purify DNA after adapter ligation. Higher yield will be obtained with high quality DNA template.
- *TransStart[®] FastPfu* Fly DNA polymerase cannot incorporate dUTP. dUTP-containing primers or templates in the reaction are not recommended.

For research use only, not for clinical diagnosis.

Version number: V5-202404

Service telephone +86-10-57815020

Service email complaints@transgen.com

