

## ProteinRuler® IV (30-200 kDa)

Cat.No. DR401

Storage: at -20°C for two years

Concentration

2 µg/5 µl for 30 kDa band; 1 µg/5 µl for 100 kDa band; 0.5 µg/5 µl for each of other bands

Description

ProteinRuler® IV is composed of six unstained recombinant proteins (40 kDa, 60 kDa, 80 kDa, 100 kDa, 150 kDa, 200 kDa) and one blue prestained recombinant protein (30 kDa). The prestained band allows monitoring electrophoresis. The 100 kDa band has doubled intensity to serve as a reference band.

- MW range from 30 to 200 kDa.
- Ready-to-use format, direct load on gels without heating.

Storage buffer

50 mM Tris-HCl (pH 6.8), 5 mM EDTA, 10 mM DTT, 10% glycerol, 1% SDS, 0.01% phenol red


Notes

- Mix well before use.
- Use 5 µl/well for mini gel, 10 µl/well for larger gel.

Electrophoresis condition

For BioRad Mini Electrophoresis installation, electrophoresis at 200 V, 50 minutes.

ProteinRuler® IV  
(30-200 kDa)



8% Tris-glycine SDS gel (5 µl/well)

**For research use only, not for clinical diagnosis.**

Service telephone +86-10-57815020

Service email complaints@transgen.com

