

EasySee® Western Marker (25-90 kDa)

Cat. No. DM201

Storage: at -20°C for two years

Description

EasySee® Western Marker is composed of eight proteins ranging from 25 to 90 kDa. The 25 kDa, 45 kDa and 85 kDa bands are prestained allowing easy identification and monitoring electrophoresis and membrane transfer. Other five bands contain several IgG binding sites, allowing marker visualization using the same reagents and protocol for your target proteins. These no-dye-attached proteins provide more accurate molecular weight estimation.

Highlights

- Three prestained bands for monitoring electrophoresis and membrane transfer.
- No label, no dye attached to other five recombinant protein bands.
- Five recombinant protein bands contain IgG binding sites, which can be developed with the standard Western blot substrates.
- Ready-to-use format, direct load on gels without heating.

Storage buffer

125 mM Tris-HCl (pH 6.8), 5 mM EDTA, 10 mM DTT, 17.4% glycerol, 3% SDS, 0.01% phenol red


Notes

- Mix well before use.
- Use 5 µl/well for mini gel, 10 µl/well for larger gel.

Electrophoresis condition

For BioRad Mini Electrophoresis installation and transfer system, electrophoresis at 200 V for 50 minutes, transfer to membrane at 200 mA for 2-3 hours.

EasySee® Western Marker (25-90 kDa)



For research use only, not for clinical diagnosis.

Service telephone +86-10-57815020

Service email complaints@transgen.com

