

# TransStart® FastPfu DNA Polymerase

Please read the datasheet carefully prior to use

Cat. No. AP221

Store at: -20°C for two years

Concentration: 2.5 units/μl

## Description:

TransStart® FastPfu DNA Polymerase is a hot start high-fidelity DNA polymerase used for fast PCR. Offering a high amplification efficiency and high extension rate (4 kb/min, it is 8 times that of ordinary Pfu enzymes), TransStart® FastPfu DNA Polymerase amends the defects of low amplification efficiency, low yield and low extension rate (0.5 kb/min) in common Pfu polymerase, and greatly shortens reaction time.

- Offers 54-fold fidelity as compared to *EasyTaq*® DNA Polymerase.
- PCR products can be directly cloned into *pEASY*®-Blunt vectors.
- Amplification of genomic DNA fragment up to 15 kb.
- Amplification of plasmid DNA fragment up to 20 kb.

## Features:

- Hot start, high specificity.
- High amplification efficiency.
- Fast and high fidelity.

## Applications

- Amplifies complex and high GC/ AT templates.
- High fidelity and fast PCR, blunt end cloning, site-directed mutagenesis.
- Amplifies long fragment.

## Kit Contents

TransStart® FastPfu DNA	AP221-01/11	AP221-02/12	AP221-03/13
Polymerase	250 U×1	500 U×1	500 U×6
5×TransStart® FastPfu Buffer	1.2 ml×1	1.2 ml×2	1.2 ml×12
2.5 mM dNTPs	- / 500 μl×1	- / 1 ml×1	- / 1 ml×6
6×DNA Loading Buffer	500 μl×1	1 ml×1	1 ml×2
50 mM MgSO <sub>4</sub>	200 μl×1	400 μl×1	1 ml×1
Complimentary Component	200 μl×1	400 μl×1	1 ml×1
PCR Stimulant			

## Storage Buffer

50 mM Tris-HCl (pH 8.2), 0.1 mM EDTA, 1 mM DTT, Stabilizers, 50% glycerol

## 5×TransStart® FastPfu Buffer with Mg<sup>2+</sup>

100 mM Tris-SO<sub>4</sub> (pH 9.2), 200 mM KCl, 10 mM MgSO<sub>4</sub>, 10% Glycerol, others

## Reaction Components (50 μl reaction volumes)



Component	Volume	Final Concentration
Template	Variable	As required
Forward Primer (10 $\mu$ M)	1 $\mu$ l	0.2 $\mu$ M
Reverse Primer (10 $\mu$ M)	1 $\mu$ l	0.2 $\mu$ M
5 $\times$ <i>TransStart</i> <sup>®</sup> <i>FastPfu</i> Buffer	10 $\mu$ l	1 $\times$
2.5 mM dNTPs	4 $\mu$ l	0.2 mM
<i>TransStart</i> <sup>®</sup> <i>FastPfu</i> DNA Polymerase	1 $\mu$ l	2.5 units
Nuclease-free Water	Variable	-
Total volume	50 $\mu$ l	-

Suggested conditions (50  $\mu$ l reaction volumes)

Template	Input
Genomic DNA	10-500 ng
Plasmid DNA	1-30 ng
cDNA	1-2 $\mu$ l cDNA from RT reaction (50-500 ng RNA for RT reaction)

PCR

Number of cycles	Temperature	Time
1 cycle	95°C	2 min
30-35 cycles	95°C	20 sec
	Tm-5°C	20 sec
	72°C	4 kb/min
1 cycle	72°C	5 min

PCR Stimulant

PCR Stimulant is used to optimize the amplification of complex templates or high GC/ AT templates. The amplification of the Pfu series of enzymes is enhanced significantly. The concentration of the storage solution is 5 $\times$ , and the concentration of the working solution can be adjusted between 0.5 $\times$ -2.5 $\times$

Notes

- For GC-rich templates, the recommended denaturation temperature is 98°C
- To ensure high fidelity, we recommend using high quality dNTPs. dNTPs containing dUTP cannot be used.
- It is recommended to add *TransStart*<sup>®</sup> *FastPfu* DNA Polymerase to the reaction system in the last step.
- If 5 $\times$  *TransStart*<sup>®</sup> *FastPfu* Buffer has a small amount of precipitation after thawing, please heat it in a 37°C water bath and mix it for use.

FOR RESEARCH USE ONLY

