

MagicPure[®] Blood Genomic DNA Kit

Please read the manual carefully before use

Cat. No. EC101

Version No. Version 1.1

Storage: at 15°C-30°C for 18 months. Avoid freezing.

Description

MagicPure[®] Blood Genomic DNA Kit provides an easy, fast, and effective method for isolating high-quality DNA. This kit is based on the specific interaction between nucleic acids and magnetic beads in the presence of specially formulated buffer.

MagicPure[®] Blood Genomic DNA Kit is designed for the purification of genomic DNA from (50 µl -250 µl) of fresh, frozen, or anticoagulant-treated blood. The extracted genomic DNA is suitable for downstream applications such as enzymatic digestion, PCR, and Southern blotting. This kit is compatible with magnetic rod-based high-throughput nucleic acid extractors.

- Simple and fast, no centrifugation required.
- High yield and purity.

Kit Contents

Component	EC101-01/11 (50 rxns)
Binding Buffer 17 (BB17)	18 ml
Clean Buffer 17 (CB17)	50 ml
Wash Buffer 17 (WB17)	12 ml
Elution Buffer (EB)	10 ml
Proteinase K (20 mg/ml)	1 ml
Magnetic Blood Beads	800 µl
Magnetic Stand (16 hole)	1 each/-

Starting material

- Whole blood, short term storage: 2-8°C for up to 1 week; long term storage: -70°C.
- Avoid repeated freezing and thawing of the whole blood (no more than three times) .

Procedures

Before starting, add the below indicated volume of 100% ethanol into the CB17 and WB17.

Component	EC101
Clean Buffer 17 (CB17)	50 ml
Wash Buffer 17 (WB17)	48 ml

All magnetic separation steps are carried out at room temperature. Mix the magnetic beads well by vortexing before use.

1. Add 50-250 µl of anticoagulated whole blood sample to a 1.5 ml microcentrifuge tube.
2. Add 300 µl of BB17 and 20 µl of Proteinase K into the microcentrifuge tube. Mix well by vortexing.
3. Incubate at 56°C for 10 minutes, and vortex 1-2 times during incubation.
4. Add 400 µl of isopropanol to the microcentrifuge tube. Mix by vortexing for 10 seconds. Add 15 µl of well-mixed Magnetic Beads into the microcentrifuge tube.
5. Vortex the microcentrifuge tube for 1 minute, and then stay still for 2 minutes.
6. Repeat Step 5 three times.



7. Place the microcentrifuge tubes onto the magnetic stand until the beads are pelleted against the magnet. Remove as much supernatant as possible, be careful not to remove any beads.
(Suggestions for beads separation: after placing the microcentrifuge tubes onto the magnetic stand, gently turn the tubes left and right to attach the beads to the magnet, then invert the magnetic stand 2-3 times to 'rinse' the tube cap with the supernatant. Incubate at room temperature for 30 seconds.)
8. Remove the microcentrifuge tubes from the magnetic stand, add 800 μ l of CB17 (make sure ethanol has been added) to the microcentrifuge tubes, then vortex the microcentrifuge tubes for 2 minutes. Place the microcentrifuge tubes onto the magnetic stand, and then discard the supernatant as in Step 7.
9. Repeat Step 8 one time.
10. Remove the microcentrifuge tubes from the magnetic stand, add 500 μ l of WB17 (make sure ethanol has been added) to the microcentrifuge tubes, then vortex the microcentrifuge tubes for 2 minutes. Place the microcentrifuge tubes onto the magnetic stand, and then discard the supernatant as in Step 7.
11. Repeat Step 10 one time.
12. Air-dry the uncapped beads on the magnetic stand for 8-10 minutes.
13. Remove the microcentrifuge tubes from the magnetic stand, add 80-200 μ l of EB to the microcentrifuge tubes. Mix gently by pipetting up and down several times to resuspend the beads and incubate at 65°C for 5-10 minutes. Mix gently by pipetting up and down once or twice during incubation.
14. Place the microcentrifuge tubes onto the magnetic stand, separate beads as in Step 7. Carefully transfer the supernatant into a clean 1.5 ml tube. Avoid collecting beads during the transfer, and store the purified DNA at -20°C.

Notes

- Use fresh blood sample and avoid repeated thawing and freezing.
- Mix the magnetic beads well by vortexing before use.
- Use sterile tubes and pipette tips to avoid the DNase contamination.

For research use only, not for clinical diagnosis.

Version number: V1.1-202512

Service telephone +86-10-57815020

Service email custserv@transgenbiotech.com

