

ProteinFind[®] Anti-GST Mouse Monoclonal Antibody

Cat. No. HT601

Storage: at 2-8°C for one month, at -20°C for one year

Description


ProteinFind[®] Anti-GST Mouse Monoclonal Antibody is a purified monoclonal antibody against yeast Y258 GST recombinant proteins that detects GST fusion proteins.

Suggested Dilution

- Western Blot: 1:500-10000 dilution.
- ELISA: 1:500-5000 dilution.
- IP: 1:100-500 dilution.

Notes

- Aliquot prior to use, and avoid repeated freezing and thawing.
- Diluted antibody should be used within 12 hours.



Detection of GST fusion protein expression by WB

FOR RESEARCH USE ONLY