

6×DNA Loading Buffer

Please read the datasheet carefully prior to use.

Cat. No. GH101

Storage: at room temperature for two years

Description

6×DNA Loading Buffer is used as loading buffer in nucleic acid electrophoresis. Prior to loading, add appropriate volume of 6×DNA Loading Buffer to DNA sample to make its working concentration at 1×, and then load the DNA samples into the gel for electrophoresis.

For research use only, not for clinical diagnosis

Version number: V1-202008

Service telephone +86-10-57815020

Service email complaints@transgen.com

