

TransCult[®] LB Agar Plate (Ampicillin)

Cat. No. CP111

Storage: at 2-8°C for three months

Description

This product is a ready-to-use LB agar plate medium containing 50 µg/ml ampicillin. The plate diameter is 9 cm, and the medium volume of each plate is 20 ml. It is mainly used for the cultivation of *E. coli*.

Kit Contents

Catalog Number	Specification
CP111-01	9 cm plate, 10 plates/box

Notes

- The plate is vacuum packed. If swelling or air leakage occurs, do not use it.
- Antibiotic is unstable, please use it as soon as possible before the expiration date.

For research use only, not for clinical diagnosis.

Service telephone +86-10-57815020

Service email complaints@transgen.com

