

TransCult[®] LB Agar Plate (Kanamycin)

Cat. No. CP121

Storage: at 2-8°C for three months.

Description

This product is a ready-to-use LB agar plate medium containing 50 µg/ml kanamycin. The plate diameter is 9 cm, and the medium volume of each plate is 20 ml.

Product Contents

Catalog Number	Specification
CP121-01	9 cm plate, 10 plates/box

For research use only, not for clinical diagnosis.

Service telephone +86-10-57815020

Service email complaints@transgen.com

