

# *ProteinFind*<sup>®</sup> Anti-MBP Mouse Monoclonal Antibody

Cat. No. HT701

**Storage:** at 2-8°C for one month, at -20°C for one year

## Description


*ProteinFind*<sup>®</sup> Anti-MBP Mouse Monoclonal Antibody is a purified monoclonal antibody that detects MBP fusion proteins.

## Suggested Dilution

- Western Blot: 1:500-10000 dilution.
- ELISA: 1:500-5000 dilution.
- IP: 1:100-500 dilution.

## Notes

- Aliquot prior to use, and avoid repeated freezing and thawing.
- Diluted antibody should be used within 12 hours.



Detection of MBP fusion protein expression by WB

**For research use only, not for clinical diagnosis.**

Service telephone +86-10-57815020

Service email [complaints@transgen.com](mailto:complaints@transgen.com)

