

Trypsin

Cat. No. FG301

Storage: at -20°C for 18 months

Description

Trypsin contains porcine trypsin (0.25%) , EDTA (+/-), and phenol red. It does not contain calcium and magnesium ion. Trypsin solution has been filter-sterilized and it can be used for cell dissociation.

Note

Aliquots after thawing; avoiding repeated freezing and thawing; store at 2-8°C for two weeks, or at -20°C for one year.

For research use only, not for clinical diagnosis.

Service telephone +86-10-57815020

Service email complaints@transgen.com

