

Trypsin

Cat. No. FG301

Storage: at -18°C or below for 18 months

Description

Trypsin contains porcine trypsin (0.25%) , EDTA (+/-), and phenol red. It does not contain calcium and magnesium ion. This product has been sterilized by filtration and tested for porcine parvovirus and mycoplasma. It features rapid digestion and low toxicity, making it suitable for the dissociation of tissues and monolayer cells.

Note

The optimal digestion time should be determined based on the characteristics of different tissues and cells. Prolonged digestion is not recommended, as it may affect cell viability.

Aliquots after thawing; avoiding repeated freezing and thawing; store at 2-8°C for two weeks.

For research use only, not for clinical diagnosis.

Service telephone +86-10-57815020

Service email complaints@transgen.com

