

Penicillin-Streptomycin (100×)

Cat. No. FG101

Storage: at -20°C for one year

Description

Penicillin-Streptomycin (100×) contains 10 kU/ml of penicillin and 10 mg/ml of streptomycin. The solution has been filter-sterilized. It can be used for cell culture at a final concentration of 1x.

Note

Aliquots after thawing; avoiding repeated freezing and thawing; store at 2-8°C for two weeks, or at -20°C for one year.

For research use only, not for clinical diagnosis.

Service telephone +86-10-57815020

Service email complaints@transgen.com

