

## *ProteinFind*<sup>®</sup> Anti-c-Myc Mouse Monoclonal Antibody

Cat. No. HT101

Storage: at 2-8°C for one month, at -20°C for one year

### Description


*ProteinFind*<sup>®</sup> Anti-c-Myc Mouse Monoclonal Antibody is a purified monoclonal antibody that detects recombinant proteins containing the c-Myc (EQKLISEEDL) epitope tag.

### Suggested Dilution

- Western Blot: 1:500-10000 dilution.
- ELISA: 1:500-5000 dilution.
- IF: 1:100-500 dilution.
- IP: 1:100-500 dilution.


### Notes

- Aliquot prior to use, and avoid repeated freezing and thawing.
- Diluted antibody should be used within 12 hours.



M: *EasySee*<sup>®</sup> Western Marker

Detection of c-Myc fusion protein expression in HEK-293 cell by WB and IF



FOR RESEARCH USE ONLY