
Material Safety Data Sheet



Revision Date: 1-June-2016

1. Product and Company Information

Product Name *Transsetta (DE3) Chemically Competent Cell*
Cat. No. CD801
Company TransGen Biotech Co., Ltd.
Technical Phone +86-400-898-0321
Fax +86-10-57815011
E-mail Address trans@transgen.com.cn

2. Hazards Identification

Not hazardous

Potential Health effects

Eyes May cause eye irritation with susceptible persons.
Skin May cause skin irritation in susceptible persons.
Inhalation May be harmful by inhalation.
Ingestion May be harmful if swallowed.

3. Composition

Chemical Name	CAS-No	Weight %
Glycerol	56-81-5	5-10

The product contains no substances which at their given concentration, are considered to be hazardous to health. We recommend handling all chemicals with caution.

4. First Aid Measures

General advice In the case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).
Eye contact Rinse thoroughly with plenty of water for at least 15 minutes, lifting lower and upper eyelids. Consult a physician.
Skin contact Wash off immediately with plenty of water for at least 15 minutes.
Inhalation Move to fresh air.

Ingestion Never give anything by mouth to an unconscious person. If symptoms persist, call a physician. Not an expected route of exposure.

Notes to physician Treat symptomatically.

5. Fire-Fighting Measures

Suitable extinguishing media Water spray. Carbon dioxide (CO₂). Foam. Dry chemical

Special protective equipment for firefighters Wear self-contained breathing apparatus and protective suit.

6. Accidental Release Measures

Personal precautions Use personal protective equipment.

Methods for cleaning up Soak up with inert absorbent material.

Environmental precautions Prevent further leakage or spillage if safe to do so.

7. Handling and Storage

Handling Always wear recommended Personal Protective Equipment. No special handling advice required.

Storage Keep in a dry, cool and well-ventilated place. Keep in properly labeled containers

8. Exposure Controls/Personal Protection

Exposure limits:

Chemical Name	OSHA PEL (TWA)	OSHA PEL (Ceiling)	ACGIH OEL (TWA)	ACGIH OEL (STEL)
Glycerol	10 mg/m ³	none	10 mg/m ³	none

9. Physical and Chemical Properties

General Information

Form: liquid

Color: transparent

Water solubility: soluble

Others: No information available.

10. Stability and Reactivity

Stability	Stable under normal conditions
Materials to avoid	No information available
Polymerization	None under normal processing.

11. Toxicological Information

Chemical Name	LD50 (oral,rat/mouse)	LD50 (dermal,rat/rabbit)	LC50 (inhalation,rat/mouse)
Glycerol	= 12600 mg/kg Oral	no data available	no data available

Principle Routes of Exposure/
Potential Health effects

Eyes	May cause eye irritation with susceptible persons.
Skin	May cause skin irritation in susceptible persons.
Inhalation	May be harmful if inhalation.
Ingestion	May be harmful if swallowed.
Carcinogenic effects	none
Mutagenic effects	none
Reproductive toxicity	none
Sensitization	none

12. Ecological Information

Ecotoxicity effects	No information available.
Mobility	See log Pow
Biodegradation	No information available.
Bioaccumulation	No information available.

Chemical Name	Freshwater Algae Data	Water Flea Data	Freshwater Fish Species Data	Microtox Data	log Pow
Glycerol 56-81-5		Daphnia species EC50=500 mg/L (24 h)			logPow-1.76

13. Disposal Considerations

Dispose of in accordance with local regulations.

14. Transport Information

IATA	Not dangerous goods
Proper shipping name	Not classified as dangerous in the meaning of transport regulations
Hazard class	none
Subsidiary Class	none
Packing group	none
UN-No	none

15. Regulatory Information

Chemical Name	TSCA
Glycerol	-

16. Other Information

For research use only. Not intended for human or animal diagnostic or therapeutic uses.

The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may present unknown hazards and should be used with caution. Since the Company cannot control the actual methods, volumes, or conditions of use, the Company shall not be held liable for any damages or losses resulting from the handling or from contact with the product as described herein. THE INFORMATION IN THIS MSDS DOES NOT CONSTITUTE A WARRANTY, EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PUPOSE

End of Safety Data Sheet